

PHLEBOTOMY ORDER OF DRAW



AEROBIC BLOOD CULTURE



ANAEROBIC BLOOD CULTURE



SODIUM CITRATE

NOTE: WHEN DRAWING A COAG TEST USING A BUTTERFLY; THE LINE MUST BE PRIMED BY DRAWING A DISCARD TUBE (SODIUM CITRATE ONLY).

ANY ORDER WITHIN GROUP



TRACE ELEMENT - SERUM (Clot Activator – Red Stripe on Label)



SERUM (CLOT ACTIVATOR)



SERUM SEPARATOR



SODIUM HEPARIN



LITHIUM HEPARIN - PST GEL

ANY ORDER WITHIN GROUP



TRACE ELEMENT – PLASMA (K2 EDTA 10.8 MG)



K2 EDTA



PPT PLASMA (WHITE PEARL TOP)



ESR 3.2% SODIUM CITRATE



K2 EDTA 5.4 MG



SODIUM FLUORIDE/ POTASSIUM OXALATE



ACD SOLUTION A